LuxaPrint Ortho Comfort

DMG Digital Enterprises SE
 Elbgaustraße 248 · 22547 Hamburg · Germany · www.dmg-dental.com
 092878/2024-04

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Produktbeschreibung

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IBOMA, alip. Urethanacrylat, alip. Urethanmethacrylat, EBPADMA, TPO, Additive

Instructions for use

English

Product description

Light-curing resin based on (meth)acrylate for the production of dental devices. For DLP/SLA printers that work at a wavelength of 385 nm.

Devices with material parameters already defined in their systems are recommended. You will find a list of validated devices on our website. Select material variants according to the wavelength used by the printer (see »Packaging«).

Intended purpose

Manufacturing of occlusal splints using 3D printing.

Limitations of use

Only insert LuxaPrint Ortho Comfort intraorally in a fully polymerized state.

Indication

Bruxism and craniomandibular dysfunction (CMD)

Contraindication

Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

People to be treated in the course of a dental procedure.

Intended user

Dentist, dental technician

Notes for use

- During handling (e.g. transferring), the material should be exposed to daylight/ambient light only for as short a time as possible, since lengthy exposure to light can lead to unwanted curing of the material.
- Wear nitrile gloves, protective clothing and goggles when working with the material/printed object.

For optimum properties of the end product:

- Ensure that the devices used during the production process are in proper condition and comply with machine and material parameters.
- Ensure clean working procedures. Residues on the devices can cause defects on the printed object.
- Use a cleaning unit with material parameters defined in the system. You will find a list of validated devices on our website.
- Allow objects to cure completely.
- Use a light-curing unit with material parameters defined in the system. You will find a list of validated devices on our website.
- Maximum curing depth during post-curing: For solid objects and light-curing on both sides, the material can be up to 7 mm thick (at a curing depth of 3.5 mm).
- Thinner layers increase the accuracy of the printed object and increase printing time.
- Position the printed object as flat as possible on the printer's build platform. A flatter positioning increases accuracy and reduces printing time.

Recommended use

Note: Observe the machine and material parameters.

Note: Maintain a processing temperature of 23 °C ± 2 °C (73.4 °F ± 3.6 °F).

Note: Scan the RFID code of the material for a higher level of process security.

Finishing the printed object

Production steps	
Recommended drip-off time after raising the platform: 10 min	
Note: Finishing should take place immediately after the construction process where possible.	
1.	Carefully remove the printed object from the build platform.
2.	Cleaning: Clean the printed object with isopropanol (approx. 99 %) using an appropriate cleaning device and specified program.
2. (Alternatively)	1. Pre-clean the printed object with isopropyl alcohol (approx. 99 %) for max. 3 min in an ultrasonic cleaner (with a brush if necessary). 2. Optional: Clean openings, boreholes and gap areas using pressurized air. Clean printed object in a separate container using clean isopropanol (approx. 99 %) for max. 2 min.
3.	Post-curing: Post-cure the printed object using the appropriate light-curing unit and the specified program.
4.	After removing the support, finish the surface of the printed object using the appropriate instruments and polishing paste.

Interactions

Residues of methacrylate-based materials (uncured material residues) may influence the setting behavior of silicone impression materials.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- Avoid skin contact with the liquid material and components prior to post-curing! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Avoid breathing in the fumes.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.
- Do not use heat-based disinfection or sterilization methods. Otherwise, the workpiece may be deformed.

Product performance characteristics

Cured material	
Shore-D hardness	≤ 72

Uncured material	
Dynamic viscosity at 50 Pa shear stress	≤ 3.0 Pa s
Light sensitivity to ambient light	≤ 5 min

Storage/disposal

- Store in a dry place at room temperature (15 – 25 °C/59 – 77 °F) and protected from light!
- Even low exposure to light can trigger polymerization.
- Always keep the container tightly sealed; carefully seal immediately after every use.
- Do not use after the expiration date!
- The cured material can be discarded together with household waste following consultation with the disposal company and in compliance with regulations.
- Incinerate the resin mixture as hazardous waste in accordance with official regulations.

Composition

IBOMA, aliphatic urethane acrylate, aliphatic urethane methacrylate, EBPADMA, TPO, additives

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